

Appendix No. 2 The Technical Specifications Related to Tender No. 2/2025/YEM

NO.	Description
ICU-1	General ICU
1	Lung Ventilator for Adult & Paediatric (with compressor or turbine)
	Ventilation Modes: volume and pressure type of breaths , Patient types Adult &paediatric .
	A/C, SIMV, CPAP/PEEP, PSV, BIPASIC , VG and Bi level ventilation modes and newer modes .
	Pressure Limited Ventilation in volume control Tadial volume 5 - 1500 ml or more than 1500 ml .
	Automatic regulation of flow in volume controlled breaths
	Non-Invasive Ventilation
	Apnea back-up ventilation active in spontaneous breathing modes.
	Flow acceleration control
	Automatic leakage compensation ensuring accurate user defined setting to the patient.
	Adjustable tube resestant compensation
	Manual inspiration
	Integrated pneumatic nebulizer facility.
	Oxygen enrichment 100% for suctioning procedure
	Measurement of resistance and compliance
	Monitoring of patient data
	Flow acceleration
	Monitoring of inspiratory O2 concentration
	Accurate measurement of flow and volume :
	Airway pressure 0- 99 cm H2O approx.
	Breathing rate:0-150 bpm or more than 150 bpm
	I:E ratio (variable) : 1:0.5 to 1:6 or better
	Flow 3 - ≥ 160 l/min or more than 160 l/min
	Peak inspiratory pressure : 0 - ≥ 120 cm H2O
	Inspiratory time 0 - ≥ 20 sec.
	Measured oxygen concentration 21-100%
	PEEP/CPAP and P low : 0 to 50 cmH2O
	Alarms :
	Power failure
	Audible and visual alarms should activate when their specific monitored parameter reaches and
	remains at the set limit to activate the alarm.
	Audible alarms should be distinct and easily identified.
	If the alarm volume is adjustable, it should not be possible to turn the volume down so low that it is not
	likely to be heard.
	Although an audible-alarm silence is acceptable, the alarm must recur automatically if the condition is
	not corrected.
	If an alarm is silenced, a visual display should clearly indicate which alarm is disabled.
	Visual alarms should be easy to identify.
	The visual alarm must be specific to the problem and remain on until the alarm condition is corrected; it
	should not be possible to turn off the visual alarm.
	All alarms should be fully explained in the operator's manual.
	The delivered oxygen or oxygen/air mixture shall be monitored with an oxygen analyzer that includes an
	alarm for concentrations outside of acceptable ranges.

	Better The analyzer shall be included within the ventilator (Pramagntic O2 cell), but if O2 sensor
	separately and placed in line with the breathing circuit, shell be supplied with 3 - 4 spare O2 sensors
	(Gelvanic O2 cell 7000 Hr)
	The delivered oxygen or oxygen/air mixture shall be heated and humidified with an add-on-heated
	humidifier.
	Time cycled, volume constant, pressure controlled ventilator >
	Integrated graphical display , at least 14 " colour monitor , for settings, values , pressure and flow
	curves.
	Internal battery operation at least 3 hr or more
	Should have trending facility for facility for 1 week or more .
	Automatic gas switch-over if oxygen pressure decreases
	The machin can be oprated by both central gase supplay and cylinder
	Digital port interface RS 232
	Standard accessories :
	Ventilator mobile trolley
	Support arm for patient hoses
	Heated humidifier complete with accessories
	Reusable autoclavable silicone breathing sets for adult &other one pediatric incluing water trap bacterial
	filter, and disposable autoclavable breathing sets (10 pcs for for adult &others for pediatric incluing
	water trap bacterial filter)
	Test lung for adult &pediatric applications
	Type of Neubilizer: built in
	1 temperature probe 1 Kit of adults hoses for humidifir
	1 O2 cylinder and medical pressure pre-set regulator.
	Operation manual &service manual
	Power requirement: 220 V/50 / 60 HZ.
	Shall comply to CE marked or FDA Certificate.
	Installation &Commissioning :
	2 Veers warrants with an ear north from the data of installation and commissioning if the CC marked or 2
	3 Years warranty with spear parts from the date of installation and commissioning if the CE marked or 2
	Years waranty with spear parts from the date of installation and commissioning if FDA Certificate
	Technical Data Sheet must accompany the offers
	User /Nurses training, by Specialist from the Supplier.
	Note:
	All equipment needing consumables must allow the possibility to use generic and/or locally made
	consumables and/or disposables. Compliance to this condition must be declared here by the bidders.
	The system offered shall be designed to operate normally under the conditions of the purchaser's
	country. The conditions include Power Supply, Climate, Temperature, Humidity, etc.
OP	Operation Threator Rooms
OP-1	Electric Cautery
	Diathermy apparatus, generator high frequency
	bipolar and monopolar channels
	Up to 350 watts
	Cutting and coagulation mode
	Operating frequency up to 350 KHz
	Microcontroller-based multi-programmable by the user
	Power supply 220-240V 50/60 Hz
	Resettable over current breaker shall be fitted for protection
	Alarm audiovisual in case of power supply disconnection.
	Power disconnection in case of loose patient plate
	Mounted on trolly, stainless steel with swivel castor wheels with brake
	Complete with the following accessories:
	Two button handle autoclavable 1
	Two button handle cold sterilization 1
	cable for plate 1
	Reusable patient plate 2
<u></u>	I resease patient plate 2

	Waterproof double foot switch 1
	Short, general surgery blade electrode
	Long, general surgery blade electrode
	Short, general surgery needle electrode
	Short, general surgery ball electrode
	Diathermy Quiver Insulated, diathermy Cable Red rubber 5 meter
	Supply with :
	User manual in English
	Service manual in English
	Shall comply to CE marked or FDA Certificate.
	Installation & Commissioning :
	2 Years waranty with spear parts from the date of installation and commissioning .
	User /Nurses training, by Specialist from the Suplier.
	Note:
	Note:
	All equipment needing consumables must allow the possibility to use generic and/or locally made consumables and/or disposables. Compliance to this condition must be declared here by the bidders. The system offered shall be designed to operate normally under the conditions of the purchaser's country. The conditions include Power Supply Climate. Temperature, Hymidity, etc.
00.0	country. The conditions include Power Supply, Climate, Temperature, Humidity, etc.
OP-2	Autoclave Sterilization Machine
	STEAM STERILIZER Vertical
	Chamber volume: 80-100 liters
	(Vacuum – Steam pressure –dry) Phases
	Sterilization temperatures : 121°C / 134°C
	Microprocessor control
	Seamless stainless steel chamber
	Over pressure Safety valve
	100- 240 V - 50 Hz
	Supplied with:
	Set of two door gasket
	set of four bacteriologoical filters
	With suitable water distiller (glass type not acceptable)
	The product should be at least CE/FDA certified high quality.
	2 Years waranty with spear parts from the date of installation and commissioning
	Technical Data Sheet must accompany the offers User /Nurses training, by Specialist from the Suplier.
OP-3	Dry Heat Sterilization Machine
	Temp.range: RT+10-300 °C
	Temp.Rsolaution Ratio:0.1°C
	Temp.Motion:± 1°C
	Temp.Uniformity:± 3.5°C
	Inner Chamber: Mirror Polished Aluminum
	Outer shell : cold rollaing steel electrostatic
	Volum:85L
	Load per rack 15Kg
	selif space 35mm
	sheilf: 2
	POWER RATING 50\60hz:690*695*985
	Inner chamber size : 450*420*450 Extrnal size: 590*610*830
OP-4	Hydraulic Operating Table for General Surgery
	T shaped base,5 sections (divided leg)
	together, easily moving wheels with brakes, base and column are constructed from st/st,
	top plate are X- ray translucent along the full length of the table, with mattress
	antistatic, head rest can be replaced
	* Working load: not less than 180 KG
	* Adjustments :
	- Height : 750 -1000 mm hydrulically using foot pedal or better
	1.700 1000 mm Hydramodny doing foot poddi of bottor

	Trandlanhura 1.1/ OF dagree manually
-	- Trendlenburg : +/- 25 degree manually
	- Back Section : -15 to + 70 degree manually
	- Side Tilting : +/- 15 degree manually
	- Leg Section : -90 degree manually
	- Head Rest : -40 to + 30 degree manually
	Castors: twin antistatic, two with foot brakes
	Accessories with necessary clamps:
	All ball joints, clamps, rods and fixations should be made from stainless steel
	Flat amount to allow to the land of the North and the North and the second discount of the
	. Flat, square head rest, (table width), infinite slope adjustment, integral foam padding, 1 piece each
	. Leg plates, split type, including integral foam padding, 1 piece each
	. Patient restraint / body strap, 1 piece
	. Arm supports, 2 pieces
	. IV rod, 1 set, adjustable height Stainless Steel.
	. Shoulder rest / support with pad, 2 pieces
	. Narrow arm table with pad, 1 piece
	."D"shaped ankle pillow, 1 piece
	. Accessory / setting clamps, circular socket, 4 pieces
	. Accessory / setting clamps, rotary, 4 pieces
	1 Anaesthesia screen with radial support
	2 Leg section supports with support
	1 Body strap
	1 Head support , detachable
	1 Set of anti-static ,divided 6 parts
	- Paediatric stirrups
	-Arm and hand table
	- Stainless steel bowl with drain for gynaecologic
OP-5	Mobile Operation Light
	Surgical Operating Mobile Lamp Led's type. Shall be very stable on 4/5 castors with antistatic rubber
	bands.
	Shadow less cold light. Microprocessor Controlled
	Adjustable lighting focus by a central handle (autoclavable)
	Controllable intensity
	Lamp head diameter shall be 200 to 500mm
	Requested luminous intensity: > 90.000 - 110.000 Lux at 1.0m
	Shall run on mains and battery . When used on Battery (DC) the decreasing of intensity shall be less
	than 20%
	Lighting autonomy shall not be less than 2 Hrs
	Lamp rod shall be adjustable in height. As an indication:
	Battery shall be equipped with an automatic charger with uninterruptible connection to mains
OP-6	Laryngoscope for Adult & Child
01-0	One Handle With Four Blades sizes 1, 2, 3, 4
	Supplied With Four Spare Bulbs .
OP-7	Anaesthetist (OT) Trolley with Drawer
<u> </u>	Fully made from Stainless Steel / Plastic or aluminium with appropriate reinforcements. Mobile on 4
	antistatic castors , 2 of which shall have foot brakes
	The unit shall carry 2 height adjustable trays, 4 drawers with silent sliders and smooth edges as well as
	locks. Offer must specify maximum load on drawers.
	Upper Part shall carry Stainless Steel / plastic or aluminium pigeon holes on all width and welded in
	elevation on specific frame. All around Stainless Steel rails with 3 hooks
	Shall be protected with efficient rubber bumpers
OP-8	Patient Monitor
	Multi-parameter monitor with 12 inch TFT displaying 6 - 8 waveforms :
	72 Hr (or more) Trends graphical and Tabular.
	Built in battery 3 Hrs (Battery is accessable).
	Arythmia analysis.

	HD/DD ECC/DEED NIDD CDO2 2TEMD Coveribble Fi CO2 Et CO2 Module
	HR/PR, ECG/RESP-NIBP-SPO2-2TEMP &available Fi CO2, Et CO2 Module
	rapid and easy to all functions and parametr settings, alarm mute, on/off freeze, on/off recordr
	HR/PR, ECG/RESP-NIBP-SPO2-2TEMP-Respiration and Temprature, alarm recall for last
	12 alarms with date and time. Built-in battery capable for 30 minutes operation - Display of Systolic, Diastolic and mean pressure
	• ECG:
	- 5/12 Lead, ECG Monitoring, I, II, III, V1,V2, V3,V4,V5,V6, AVR ,AVF, AVL
	• Respiration:
	- Display of respiration waveform.
	- Apnea Alarm.
	• SP 02
	- Range: 0-100%.
	- Pulse Rate range: 25-250 bpm
INST	Surgical Instrument
INST-1	Delivery Set
11401-1	Each Set Consists of :
	Bowl with sterile swabs 2 pc
	Sponge holder 2 pc
	Urinary catheters 2 pc
	Arteries forceps(Curved, Straight) 2 pc
	Kocher Forceps (Curved, Straight) 2 pc
	Episiotomy Scissors 2 pc
	Dissecting Forceps(Tooth , Plain) 2 pc
	Needle Holder 2 pc
	Cord cutting 2 pc
	Kidney Tray 2 pc
	Surgical suture 1 pc
	Strile draping sheet 2 pc
	Sterile gloves 2 pc
	Swabs/gauze 2 pc
	Gynae pads 2 pc
	Syringe(10cc) 2 pc
	Needles 2 pc
INST-2	General Surgical Set (Large)
	Each Set Consists of:
	Bowl Lotion 250 mm 2 pcs
	Gallipot/Cup S/S 180 ml 2 pcs
	Gallipot/Cup S/S 250 ml 2 pcs
	Tray Kidney shape with Lid S/S 200 X 30 mm 2 pcs
	Tray Kidney shape with Lid S/S 300 X 50 mm 2 pcs
	Forceps artery Kocher COF BJ.S/S 180mm. 20 pcs
	Forceps Artery Spencer Wells CO Flat (bj).S/S 200 mm 15 pcs
	Forceps Dissecting Plain serrated jaws S/S 150 mm 6 pcs
	Forceps Dissecting Plain serrated jaws S/S 250 mm 1 pc
	Forceps Dissecting narrow, serrated jaws S/S 200 mm 1 pc
	Forceps Dissecting Lane's 2x3 teeth 180 mm 3 pcs
	Forceps Dissecting Waugh, fine serr.jaws 1x2 Tth 200 mm 1 pc
	Forceps Tissue Allis fine 4 x 5 Teeth BJ.S/S 150 mm 4 pcs
	Forceps Towel Holding cross Action S/S 100mm 8 pcs
	Forceps Towel Moynihan, Double Teeth (bj) 190mm 4 pcs
	Forceps Sponge Holding Rampley BJ S/S 240 mm 5 pcs
	Handle size 3 for Scalple blades size 10-15 1 pc
	Handle size 4 for Scalple blades size 20-24 1 pc
	Needle Holder Kilner's 165 mm 1 pc
	Needle Holder Mayo Hegar BJ.S/S 180mm 2 pcs
	Scissors Blunt Point Straight SJ. S/S 180mm 1 pc
	Scissors Mayo's CO Flat 200mm S/S 1 pc
_	

_	lo : N
	Scissors Mayo's Straight 200mm S/S 3 pcs
	Scissors McIndoe's CO Flat 190mm S/S 1 pc
	Retr.Deaver #1-Blade 25mm, Length 305mm 1 pc
	Retr.Deaver #5-Blade 75mm, Length 305mm 1 pc
	Retr. Lange beck, H/H 26x64 Blade, 210 mm Long 2 pcs
	Retr. Morris Double-ended (52x38)&(52x64) mm 2 pcs
	Retr. Walton, malleable copper 13mm length 305mm 1 pc
	Needle Aneurysm S/S . Small 146 mm 1 pc
	Pin Safety Mayo, forceps holding, 120mm 4 pcs
	Suction Tube McIndoe 1 pc
	Suction Tube Poole's S/S 225mm 1 pc
	Sterilizing and storing case St. St. with filters 260x150x50 mm 2 pcs
	Wire Cutting, Tungsten .Carbide inserts 180mm. 1 pc
INST-3	Instrument Trolley with Stainless Steel Shelves
	Fully made from Stainless Steel .
	Mobile on 4 antistatic castors, 2 of which shall have foot brakes
	The unit shall carry 2 height adjustable trays with rebatted edges
	Indicative Dimensions: 900 x 600mm
MNCH	MNCH Departemnt
IVINCH	-
MNICH	Ultrasound machine\Digital Ultrasonic Diagnostic Imaging System\Powerful B/W System with
MNCH-1	Full Screen; Exceptional Ergonomic Design ; Enhanced Diagnostic Con dence and Streamlined
	Work + Ultrasound Printer + UPS 2kVA
	Digital Ultrasonic Diagnostic Imaging System
	1. General Description
	The offered equipment shall be a high-quality Digital Ultrasonic Diagnostic Imaging System (B/W)
	suitable for use in a variety of clinical settings. The system must provide excellent image quality,
	streamlined workflow, ergonomic design, and reliable performance for daily diagnostic use.
	Technical Specifications
	Digital Ultrasonic Diagnostic Imaging System (B/W)
	Full-screen display with high-resolution imaging
	2D, M-Mode, Zoom, Measurement &Calculation Functions
	Adjustable monitor, control panel, and user-friendly interface
	Internal storage with USB export capability
	Supports multiple probes (Convex, Linear, Transvaginal – as per clinical need)
	Enhanced image processing and clarity for improved diagnostic accuracy
	Fast boot-up, easy patient data management, and quick image access
	Compatible with UPS for uninterrupted operation
	Standard Accessories Included
	High-resolution thermal printer for direct image printing
	Uninterruptible Power Supply (Minimum 2 kVA capacity)
	One Convex Probe and Linear prob (standard); others optional as required
	Complete operation and service manual in English
	Included for standard installation and operation
MNCH-2	Fetal Doppler
	Multi-purpose pocket Doppler for use within general practices and health centres . The Unit shall be
	supplied with plug-in version transducer to provide a versatile diagnostic aid for foetal heart rate (display
	FHR)
	Unit shall have only two controls; an on/off switch on the transducer and a volume (1 W Audible Out
	Put) control on the main unit
	Transducer: 2 MHz - 2.5 MHz ± 10% Waterproof.
	Power shall be from either a conventional 9 volts battery or an rechargeable battery
	Battery indicator shall provide advance warning of low power levels
	Shall be supplied in a compact carrying case fitted with a shaped interior to securely hold equipment
	and accessories
-	Shall be supplied with pocket / belt clip for convenient location during or between procedures.
-	Supplied with:
-	. 1 transducer 2 MHz ± 10%
	. (I'alisuuce)

	. 1 carrying case
	. Gel, electroconductive. 2x250g
	. 2 rechargeable batteries (if the model proposed functions with rechargeable batteries)
	Battery charger (if the model proposed functions with rechargeable batteries)
MNCH-3	Patient Wheelchair
	Metal construction, Nicle chrome plated, foldable, with two wheels 25-inch dia, and two 8 inch
	dia, and hand brakes, aluminium foot rest, seat and back are blue leather, adult size. Offer must specify
	maximum load on chair.
MNCH-4	
	Weighing Infant Scale vveigning range. up เอ zo kg.
	Minimum graduation: 5 g.
	Accuracy: ±5 g.
	Precision: ±5 g.
	Display: kg.
	Tare function.
	Autohold function.
	Automatic switch-off.
	Auto-calibration with each switch-on.
	Large LCD display.
	Reading time max 5 seconds.
	Splash proof and shock resistant.
	Smooth surface/finishing for easy cleaning/disinfection. Material: Body and tray made of non absorbent ABS plastic.
	Colour: white or light grey.
	With or without levelling feet.
MMCH-5	Infant Warmer
IVIIVICIT-3	Infant warmer to be used in neonatology.
	Infant warner to be used in heariatology.
	Microprocessor control system and manual controlling the temp between skin and air automatically .
	Specifications
1	OPECHICALIONS
	, •
	· Temperature control:
	· Temperature control: · Range 30-38° C
	· Temperature control:
	Temperature control: Range 30-38° C Skin range 25 – 42 ° Increment 0.1°
	· Temperature control: · Range 30-38° C · Skin range 25 – 42 °
	Temperature control: Range 30-38° C Skin range 25 – 42 ° Increment 0.1° Monitor display for all parameters
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.)
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual)
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect
	• Temperature control: • Range 30-38° C • Skin range 25 – 42° • Increment 0.1° • Monitor display for all parameters • Control Unit (to be supplied with.) • Automatic heat control type • Set point mechanism • Heater Indicator. • built-in halogen or led light • APGAR score facility • Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) • High air temperature • Sensor disconnect • Power Failure
	• Temperature control: • Range 30-38° C • Skin range 25 – 42° • Increment 0.1° • Monitor display for all parameters • Control Unit (to be supplied with.) • Automatic heat control type • Set point mechanism • Heater Indicator. • built-in halogen or led light • APGAR score facility • Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) • High air temperature • Sensor disconnect • Power Failure • System fail.
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes:
	• Temperature control: • Range 30-38° C • Skin range 25 – 42 ° • Increment 0.1° • Monitor display for all parameters • Control Unit (to be supplied with.) • Automatic heat control type • Set point mechanism • Heater Indicator. • built-in halogen or led light • APGAR score facility • Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) • High air temperature • Sensor disconnect • Power Failure • System fail. • Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: • Self- check features
	• Temperature control: • Range 30-38° C • Skin range 25 – 42° • Increment 0.1° • Monitor display for all parameters • Control Unit (to be supplied with.) • Automatic heat control type • Set point mechanism • Heater Indicator. • built-in halogen or led light • APGAR score facility • Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) • High air temperature • Sensor disconnect • Power Failure • System fail. • Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: • Self- check features • Breaks for casters
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: Self- check features Breaks for casters Skin sensor
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: Self- check features Breaks for casters Skin sensor Shelf &I.V Pole
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: Self- check features Breaks for casters Skin sensor Shelf &I.V Pole Protection against breaks and bursts of radiant and light source
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: Self- check features Breaks for casters Skin sensor Shelf &I.V Pole Protection against breaks and bursts of radiant and light source Spares and all accessories
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: Self- check features Breaks for casters Skin sensor Shelf &I. V Pole Protection against breaks and bursts of radiant and light source Spares and all accessories Power Supply:
	Temperature control: Range 30-38° C Skin range 25 – 42° Increment 0.1° Monitor display for all parameters Control Unit (to be supplied with.) Automatic heat control type Set point mechanism Heater Indicator. built-in halogen or led light APGAR score facility Accessible from 4 sides with plexi glass safety frame Alarms (Audible and Visual) High air temperature Sensor disconnect Power Failure System fail. Alarm in manual mode: every 15 minutes with automatic shutoff The warmer should includes: Self- check features Breaks for casters Skin sensor Shelf &I.V Pole Protection against breaks and bursts of radiant and light source Spares and all accessories

	· Backup thermostat
	Temperature Reusable sensor (skin 1)
	Standards, Safety and Training :
	Should be FDA or CE &TUV approved product
	User /Nurses training, by Specialist from the Supplier. In the first of Communication in the Supplier. The first of Communication is a second of the Supplier. The first of Communication is a secon
	• Installation &Commissioning
	2 Years waranty with spear parts from the date of installation and commissioning
	Documentations to be provided :
	User/Technical/Maintenance manuals to be supplied in English.
LAB	Laboratory and blood room
LAB-1	Automatic Hematology Analyzer (CBC) – Complete Blood Counter
	Principle: Electrical Impedance Method- WBC / RBC / PLT, Colorimetric- HGB
	Open or Closed System, Reagents device must have an economic consumption of reagents mechanism.
	Reportable Parameters: 5 Part Differential, 23 Parameters - WBC, LYM#, MID#, GRA#, LYM
	MID%,GRA%, RBC,HGB, HCT,MCV, MCH, MCHC, RDW-SD, RDW-CV, PL MPV, PDW, PCT,P-LCR,
	WBC Histogram, RBC Histogram, PLT HistogramT
	Operating Platform: Windows interface with mouse control (Not Compulsory).
	Speed: 40-60 samples per hour
	Calibration: Both Automatic and Manual Calibration mode. 3 sets of Calibration co -efficient for Whole
	Blood, Anti-coagulant Peripheral Blood and Pre-diluted Peripheral Blood.
	Quality Control: 3 QC Method; and X-R QC with appropriate QC Charts, Westgard
	Results: Internal memory can store at liest 40,000 Test values with Histograms
	Printer / Display: Builtin thermal printer with optional external printer, Color LCD, 9"
	Interface: 2- 4 USB, 1 RS232 serial port, 1 Lan, 1 PS2 Mouse, 1 PS2
	Power Supply / Input Power: 220 / 50 Hz
	Backup suitable UPS
	Supplied with:
	Operation Manual
	Service Manual
	Service/ maintenance and technical documents inclusive of schematics, component diagrams, trouble
	shooting, spare parts ordering information etc.
	2 Years comprehensive warranty, from the date of installation and commissioning
	Installation &Commissioning must be done by manufacturer engineer
	Users training, by Specialist from the Suplier.
	Shall comply with CE marked or FDA Certificate.
	NOTE
	Machine with close Or open system, the supplier must supply free quantity of reagents with machine and must mention how many free kits and tests of reagents would be supplied. The cost of soluations
-	&consumbles per Test for one Year should also be offered.
	Machine with close Or open system, should be put cost of soluations, consumbles per Test ,,should be
	put cost of reagents kits for at least 5 years (price / kit) and how many tests per kit. Device must be
	supplied with manufacturer reference test sample for calibration and quality assurance.
	Depended of above we well choise the system The system effected shall be designed to energia permally under the conditions of the purchaser's
	The system offered shall be designed to operate normally under the conditions of the purchaser's
	country. The conditions include Power Supply, Climate, Temperature, Humidity, etc.
LAB-2	Fully Automatic Chemistry Analyzer
ļ	Fully Automatic Clinical Chemistry Analyzer
	1. General Description
	Compact, benchtop, fully automatic discrete chemistry analyzer.
	2. System Performance
	Throughput: Minimum 50 -100 tests/hour.
	Analysis Types: End-point, kinetic, two-point, and optional ISE.
	Reaction Volume: 100–300 μL (programmable).
	Sample Volume: 2–50 μL (programmable).
	Reagent Positions: Minimum 20 on-board reagent slots.

	Sample Positions: Minimum 20 sample positions with continuous or batch loading.
	3. Analytical Features
	Photometric system with LED or halogen light source.
	Wavelength range: 340–800 nm (minimum 8 wavelengths).
	Temperature-controlled incubation at 37°C ± 0.2°C.
	Disposable or washable cuvettes.
	Supports full biochemistry panels including liver, renal, lipid, glucose, enzyme, and electrolyte tests.
	4. Sample and Reagent Handling
	Barcode system for sample and reagent identification (preferred).
	Clot detection and liquid level sensing.
	STAT sample access and emergency testing.
	5. Calibration and Quality Control
	Automatic multi-point calibration.
	Integrated QC with Levey-Jennings and Westgard rules.
	Long-term data storage for calibration and QC history.
	6. User Interface
	Touchscreen or PC-based interface.
	Multilingual user-friendly software (English required).
	Minimum storage for 10,000 patient results.
	Bi-directional LIS compatibility.
	7. Connectivity and Data Management
	USB and LAN ports for data transfer and connectivity.
	HL7 protocol support.
	Export formats: PDF, Excel, CSV.
	8. Power and Environmental Requirements
	Power Supply: 220–240 V AC, 50/60 Hz.
	Ambient Temperature: 15°C to 30°C.
	Humidity: Less than 80% non-condensing.
	Water type: Deionized or distilled water, if required.
	9. Accessories and Consumables
	Starter kit with reagents
	consumables for 500 tests Option
	Maintenance kit and spare cuvettes for 12 months.
	•
	UPS (minimum 2 kVA with 30-minute backup). 10. Certifications and Compliance
	· ·
	CE-IVD or FDA certification.
	Manufacturer must be ISO 13485 certified.
	Compliance with IEC electrical and EMC standards.
	11. Warranty and After-Sales Support
	Minimum 2-year warranty.
	On-site installation, calibration, and training included.
	Service response within 48 hours.
LAB-3	Fully Automated Benchtop CLIA Analyzer
	Fully Automatic CLIA (Chemiluminescence Immunoassay) Analyzer
	1. General Description
	Fully automatic bench-top chemiluminescence immunoassay (CLIA) analyzer.
	Performs qualitative and quantitative immunoassays using chemiluminescent detection.
	2. System Performance
	Throughput: Minimum 50 -80 tests/hour.
	Principle: Magnetic particle-based chemiluminescence.
	Assay Time: <30 minutes per test.
	Loading: Continuous or batch sample loading.
	Sample Types: Serum, plasma, whole blood .
	3. Sample and Reagent Handling
	Sample Positions: Minimum 15- 30 sample positions.
	Reagent Positions: Minimum 15 -30 reagent positions with onboard cooling.
	• • • • • • • • • • • • • • • • • • •

	Reagent recognition via barcode (preferred).
	Clot detection and liquid level sensing.
<u> </u>	4. Test Menu and Capabilities
	Must support assays for: TSH, T3, T4, FT3, FT4, LH, FSH, Prolactin, Beta-HCG, Ferritin, Vitamin D,
	Cortisol, PSA, and HbA1c.
	test menu including infectious disease markers, tumor markers, cardiac markers, etc.
	Test menu including SARS-CoV-2 IgM/IgG and hepatitis panels.
	5. Calibration and Quality Control
	Automatic calibration with multi-point curves.
	•
	Internal quality control.
	Graphical QC display with historical data storage.
	6. User Interface and Software
	Touchscreen or PC interface with intuitive operation.
	Software must support data management, patient records, and maintenance logs.
	Minimum 5,000 patient result storage.
	LIS compatibility (unidirectional/bidirectional).
	7. Connectivity and Data Management
	USB and LAN interfaces for data export and connectivity.
	Export formats: PDF, Excel, CSV.
	LIS/EMR integration with HL7 protocol preferred.
	8. Environmental and Power Requirements
	Power: 220–240 V AC, 50Hz.
	Ambient Temperature: 15–30°C.
	Humidity: <80% non-condensing.
	Onboard cooling for reagent storage (2–8°C).
	9. Accessories and Consumables
	Starter pack with reagents
	consumables for at least 300 test option.
	Printer (built-in or external).
	UPS (1.5–2 kVA) for backup power support.
	10. Certifications and Compliance
	CE-IVD, or FDA, certification.
	Manufacturer must be ISO 13485 certified.
	System must comply with electrical and EMC safety standards.
	11. Warranty and Support
	Minimum 2-year warranty including parts and labor.
	On-site training and installation.
	Remote support and guaranteed spare parts availability for at least 5 years.
LAB-4	Laboratory Microscope, Binocular
	Standard Laboratory Use Type.
	With Binocular Observation Tube, Mechanical Stage, Condenser And Filter Holder.
	Objectives Magnification: 4X, 10X, 40X &100X Achromatic Objective lines
	Eyepieces:Two eyepeices 10X widefield &2 eyepiece 16X.
	Stand Focusing Mechanism: built-on graduated mechanical stage fine adjustments: single lateral with
	coarse and coaxial adjustment graduated both sides
	Square stage Type
	Illumination: Condenser type: focusable condenser
	Low voltage Light source: Halogen or LED illumination System (Easy to replace the lamp)
	Tension Adjustment: Coaxial, X-Y controls
	Full range objectives magnification
	Complete with:
	2 pcs. Ground glass filter.
	3 pcs. spare lamps (ONLY FOR HALOGEN LAMP TYPE)
	2 pcs, 150 ml Bottle of immersion oil
	1 pc. Dust cover.
	1 Cabinetwith lock and ke
	1 Cleaning brush
	•

_	Complicativity Complete Assessment
	Supplied With Complete Accessories.
	• Guarantee Spare Parts: (Ensure the supply of spare part for a Period of (7-10) years from the supply
	date of the device)
	Power Supply:-
	- 100 to 240 V ~ ±10%, 50/60 Hz , 6V, 20W Halogen Bulb Illumination System or LED type
LAB-5	Bench Top Centrifuge
	Table Top Centrifuge.
	Maximum Revolution = 5,500 RPM.
	Full Automatic Microprocessor Control.
	Timer = From 0 - 60 Minute With Hold Mode .
	Capacity = 15ml x 16 Tubes.
	Supplied With Complete Accessories.
	Supplied With One Year Recommended Spare Parts.
	1 Year comprehensive warranty, from the date of installation and commissioning
	User / Nurses training, by Specialist from the Suplier.
	Shall comply to CE marked or FDA Certificate
EMG	General Emergency
EMG-1	Biphasic Defibrillator – Adult/Pediatric, with Heart Rate Alarm & Sync
	High resolution LCD monitor, ECG recorder and ECG monitor, separate ECG input, freeze function, hight
	rsolution printer, synchronous and asynchronous modes, internal rechargeable battery, visual and audible
	alarms
	Technical Specifications:
	Max. output : 260 - 360 joule
	Max.charging time to 360 joule : 5-8 seconds
	ECG monitoring time : 10 hr
	Battery charging time to full : 4 hr
	With Automated External Defibrillation (A.E.D)Mode with it's conections & Electrode
	System and malfunction alarms with visual and audible alarm indicators
	Supplied with Reusable adult pedals and children pedals &trolley
	Supplied with Nedsable addit pedals and children pedals attolley
	should be supplied with all accessories &all conection for adult chilled and another one for paediatrics
	Must be Certified with CE/FDA Markets high Quality.
	2 Years waranty with spear parts from the date of installation and commissioning
	2 Tears waranty with spear parts from the date of installation and commissioning
	Technical Data Sheet must accompany the offers User /Nurses training, by Specialist from the Suplier.
EMG-2	Portable Suction Pump
EIVIG-2	Vacuum range 0- 600mm Hg,
	, ,
	Flow rates – open flow > 20 litres/ minute
	Overflow safety protection
	Re-usable collection Jars dual not less than (1500cc) complete with caps and tubing (with disposable
	liners)
	Splash proof protection
	I be draw babis filters (bastavia 9 askassat) associate assistance and a satisface of Gadditional Filters and
	Hydrophobic filters (bacteria &exhaust) supplied with extra spares at least 2 additional Filters each.
	Adjustable vacuum level control knob (high vacuum)
	Suction capacity : 20 l/m
	Containers: 2 x not less than 1.5 litters each (Autoclavable)
	Hose : 2 m silicon rubber (Autoclavable)
	All accessories – i.e., Catheters, etc
EMG-3	ECG Recorder – 12 Channels with Complete Accessories
	12 Leads Simultaneous Acquisition 12 Channel Electrocardiograph.
	Recording channels: 3,6 &12 channels Recording ECG leads:12
	Shall have a LCD display.
	HR Thermal Printer.
	Digital filtering
	Selectors: I,II,III, aVR, aVL, aVF, V1,V2, V3, V4, V5, V6
	Input Imoedance: 50 MΩ Minimum
L	

	Frequency: 0.05 Hz - 100 Hz (-3dB)
	Time: 3.2 seconds
	Filter: Muscle filter at 36 Hz (-3 dB)
	Patient Safety Sensitivity 2 Francis / 10 mm/m/ 20 mm/m/ nucle button
	Sensitivity: 2.5mm/mV, 5mm/mV, 10 mm/mV, 20 mm/mV push button
	Auxiliary Input: 10mm deflection for 30 mV input
	Calibration: 1 mV push button effective at input amplifier
	Sequential lead print-out on 50mm width thermosensitive paper
	Print-out format: 3 Channel, 6 channel, 12 Channel
	Shall permit Automatic/Manual printing modes.
	Mains/Battery operated, Battery can sustain 30 - 60 mins.
	Power supply: 100-220V/50Hz, Universal Power Cord with 3 PIN, Must be supplied with external
	Earthing connection.
	The system offered shall be designed to operate normally under the conditions of the purchaser's
	country. The conditions include Power Supply, Climate, Temperature, Humidity, etc.
	Please declare in detail compliance of this item offered with any relevant quality and safety standards.
	ACCESSORIES
	Case, carrying
	4 Packs, thermosensitive. 50mmx20m
	4 pcs Electrodes, Limb. clamp type
	3 pcs Electrodes, Suction
	500 pcs Electrodes, Resting. Disposable
	10 pcs Adaptors, Clip.Alligator.4mm
	Gel, electroconductive. 5x250g
	1 Patient cable. Banana plug
	1 Spare rechargeable batteries
	Installation &Commissioning :
	2 Years comprehensive warranty, from the date of installation and commissioning
	User / Nurses training, by Specialist from the Supplier.
	Shall comply to CE marked or FDA Certificate .
	NOTE
	All equipment needing consumables must allow the possibility to use generic and/or locally made
	consumables and/or disposables. Compliance to this condition must be declared here by the bidders.
EMG-4	IV Stand – On Castors, Double Hook, Adjustable Height
	Drip Stand or Rod for I.V. procedures
	Shall be mobile on very stable stand with 4/5 castors, 4 cross-shaped / 5 stars-shaped. Castors shall be
	orientable and shall have a smooth running
	Rod made of chrome-plated steel tubing, Adjustable height
	Rod equipped with 2 hooks configuration, made of polyamide
EMG-5	Patient Stretcher – Removable Top with Side Rails and IV Stand
	stainless steel construction, painted with elecrostatic powder with removable top and mattress protected
	by Durable cover, anti-bacterial, anti-static, acid-resistant and waterproof for easy Clean
	4 castors five 5 diameter two of them with brakes
	Guiding Dimentions of trolley: 190 * 60 * 75 cm
EMG-6	Mercury Sphygmomanometer (Desk Type, 300 mm) – Adult & Pediatric Arm Bands
	Sphygmomanometer, mercury deskType
	Range of measurement shall be 0 to 300mmHg as standard with at least a precision of +/-3mmHg
	(static measure). Display of digits shall be easily readable (4mm high minimum)
	Mercury tank shall be monobloc and mercury proof. Safety tap shall be compulsory
	mentally tall so monosite and mercary proof delety tap onding to computerly
	Cuff shall be made of PVC, and shall ensure a constant and precise decompression, Velcro fixation
<u> </u>	. 1 Cuff for Adult
	. 1 Cuff for Children / Pediatrics
EMG-7	Resuscitation Set (Adult)
LIVIG-1	inconstitution out (Addit)

	Resuscitation Adult including:
	(Bag resuscitator with: patient valve, relief valve, oxygen reservoir Tube and face mask suitable for adults)
	Supplied with Spare Valve
EMG-8	Resuscitation Set (pediatric)
	Resuscitation Peditric including:
	(Bag resuscitator with: patient valve, relief valve, oxygen reservoir Tube and face mask)
	Supplied with Spare Valve
EMG-9	Oxygen Cylinder – 40 Liters
	Regulator O2 with Flowmeter and Humidifier
	Flow range: From 0.15 L/min to 15 L/min
	Maximum possible flow: 35 L/mn above range
	Input pressure: Not less than 3.5 bar +/- 0,7
	Suitable for all types of fitting: wall outlet, trolley.
	Autoclavable, unbreakable
	Clear vision of bubbles
	Humidification at 4 L/mn - 20° C: ~3 ml/h
	Humidification at 12 L/mn - 20° C: ~9 ml/h
	Capacity: 100 ml (distilled water)
	Porosity of diffuser: No more than 70 microns
	Humidifier, Autoclavable, for multipurpose uses in Medical and Surgical Sections. Shall have the capability to be associated with Flowmeters (O2) or Air , O2 Reserve Cylinders and O2 Generators
	Fixation to flowmeters or any compatible devices shall be easy and fast (1/4 twist shall be acceptable)
	CE marked
	Supplied With Tubing And Face Mask suitable for adults.
EMG-10	Finger Pulse Oximeter